

	A	B	C	D	E	F	G	H
	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	Gas-Plasma/VHP Compatibility	Additional Comments
1								
2	Laparoscopic Trays							
3		DDLDP-1123	Double Decker Lap Tray	Anodized Aluminum	Y	Y	Y	
4		CYSTO-224B	Cysto Tray	Anodized Aluminum	Y	Y	Y	
5		DDLDP-4113	Poly Double Decker Tray	Radel	Y	Y	Y	
6		SDLP-1130	Extra Long Lap Tray	Radel	Y	Y	Y	
7		51527	Double Decker Bariatric	Anodized Aluminum	Y	Y	Y	
8		31127 DV	Da Vinci Arm Sterilization Case	Anodized Aluminum	Y	Y	Y	
9		300605-DVC	Da Vinci Lens Tray	Anodized Aluminum	Y	Y	Y	
10								
11	Custom Procedure Trays							
12		30320-H	2 Level Tray System	Anodized Aluminum	Y	Y	Y	
13		31421-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
14		31123-T	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
15		31123-T1	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
16		31127-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
17								
18	Instrument Brackets for Custom Trays							
19		P-01	Post	Anodized Aluminum/Silicone	Y	Y	Y	
20		Tri-CUT BRACKET	1"universal bracket 3, 5 & 10 mm	Anodized Aluminum/Silicone	Y	Y	Y	
21		56-UB	5" 6 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
22		710-UB	7.5" 10 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
23		1015-UB	10" 15 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
24								
25	Instru-Trays							
26		4116-H		Polypropylene	Y	Y	Y	
27		4106H		Polypropylene	Y	Y	Y	
28		41013H		Polypropylene	Y	Y	Y	
29		41116H		Polypropylene	Y	Y	Y	
30		41317H		Polypropylene	Y	Y	Y	
31		41520H		Polypropylene	Y	Y	Y	
32		HCO-1309		Stainless Steel	Y	Y	Y	
33		HCO-1510		Stainless Steel	Y	Y	Y	
34		HCO-1711		Stainless Steel	Y	Y	Y	
35		HCO-1912		Stainless Steel	Y	Y	Y	
36								
37	Metal Micro Trays							
38		184014	Base & Cover	Anodized Aluminum	Y	Y	Y	
39		184016	Base & Cover	Anodized Aluminum	Y	Y	Y	
40		184019	Base & Cover	Anodized Aluminum	Y	Y	Y	
41		184021	Base & Cover	Anodized Aluminum	Y	Y	Y	
42		181453	Solid (no holes) Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
43		184453	Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
44		181463	Base,Cover & Holder	Anodized Aluminum	Y	Y	Y	
45		182500	Stainless Steel - 5 Clip	Stainless Steel	Y	Y	Y	
46		182062	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
47		182068	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
48								
49	Hinged Top Micro Trays							
50		30732-H	Taupe	Polypropylene	Y	Y	Y	
51		30842-H	Taupe	Polypropylene	Y	Y	Y	
52		31042-H	Taupe	Polypropylene	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	Gas-Plasma/VHP Compatibility	Additional Comments
53		31055-H	White	Polypropylene	Y	Y	Y	
54		31462-H	Taupe	Polypropylene	Y	Y	Y	
55		31772-H	Blue	Polypropylene	Y	Y	Y	
56		50Z900	White	Polypropylene	Y	Y	Y	
57								
58	Specialty Micro Trays							
59		30109HN	ER Hand Trays	Anodized Aluminum	Y	Y	Y	
60		30109CC	Camera Case	Anodized Aluminum	Y	Y	Y	
61		30011KE	Key Elevators Box	Anodized Aluminum	Y	Y	Y	
62		30011CU	Currettes Box	Anodized Aluminum	Y	Y	Y	
63		30011RS	Rhoton Suction Box	Anodized Aluminum	Y	Y	Y	
64		30011BS	Bipolar Set Box	Anodized Aluminum	Y	Y	Y	
65								
66	Metal Hinged Top Micro Trays							
67		31010-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
68		31010-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
69		31015-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
70		31015-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
71		31020-H	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
72								
73	Clear-Top Micro Trays							
74		0278	Base & Cover	Ultem	Y	Y	Y	
75		02781	Base & Cover	Ultem	Y	Y	Y	
76		027814	Base & Cover	Ultem	Y	Y	Y	
77		0625	Base & Cover	Ultem	Y	Y	Y	
78		062514	Base & Cover	Ultem	Y	Y	Y	
79		6540	Base & Cover	Ultem	Y	Y	Y	
80		7540	Base & Cover	Ultem	Y	Y	Y	
81		754015	Base & Cover	Ultem	Y	Y	Y	
82		100675	Single Decker Twin	Ultem	Y	Y	Y	
83		100615	Double Decker Twin	Ultem	Y	Y	Y	
84		100615B	Insert Tray for 100615	Ultem	Y	Y	Y	
85		101525	Base & Lid	Ultem	Y	Y	Y	
86		101525B	Full Insert	Ultem	Y	Y	Y	
87								
88	Silicone Plug-InSecure-It							
89		A237591K	One 4mm slot, two plugs	Silicone	Y	Y	Y	
90		A237592K	Two 4mm slots, two plugs	Silicone	Y	Y	Y	
91		A237593K	Three 4mm slots, two plugs	Silicone	Y	Y	Y	
92		A237593	Three 1mm slots, two plugs	Silicone	Y	Y	Y	
93		A237596	Six 1mm slots, two plugs	Silicone	Y	Y	Y	
94		A387596	Six 1mm slots, three plugs	Silicone	Y	Y	Y	
95		A3875912	Twelve 1mm slots, three plugs	Silicone	Y	Y	Y	
96		A387592	Two 4mm slots, three plugs	Silicone	Y	Y	Y	
97		A387593	Three 4mm slots,three plugs	Silicone	Y	Y	Y	
98		A387594	Four 4mm slots, three plugs	Silicone	Y	Y	Y	
99		A387595	Five 4mm slots, three plugs	Silicone	Y	Y	Y	
100		A2616901	One 12mm slots, two plugs	Silicone	Y	Y	Y	
101		A26216902	Two 8mm slots, two plugs	Silicone	Y	Y	Y	
102		A32759P2	One 5mm slot, three plugs	Silicone	Y	Y	Y	
103		A327593	Three 5mm slots, three plugs	Silicone	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	Gas-Plasma/VHP Compatibility	Additional Comments
104		A327593L	Two 12mm slots, three plugs	Silicone	Y	Y	Y	
105		A32759P1	One 12mm slots, three plugs	Silicone	Y	Y	Y	
106								
107	Scope Trays							
108		30126	Bariatric - Single Scope Case	Anodized Aluminum	Y	Y	Y	
109		30630	Bariatric - Multi Scope Case	Anodized Aluminum	Y	Y	Y	
110		30111-SSC	Scope Tray	Stainless Steel	Y	Y	Y	
111		30118-SSC	Universal Scope Case	Stainless Steel	Y	Y	Y	
112		31042-22	Hinged Top Arthro Tray	Polypropylene	Y	Y	Y	
113		3-4110	Clear Top Arthro Tray	Radel	Y	Y	Y	
114		SP-3014	Osteotome Tray	Anodized Aluminum	Y	Y	Y	
115		3-5120	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
116		3-5130	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
117		75.0047.9M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
118		75.0083.1M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
119		75.0047.5M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
120		750083	Wire Scope Tray	Stainless Steel	Y	Y	Y	
121		31772-22	Economy Scope Tray	Polypropylene	Y	Y	Y	
122		261615	Clear Top Scope Tray	Radel	Y	Y	Y	
123		302215	Scope Tray B/L 1 set Secur-Its	Radel	Y	Y	Y	
124		31462-22	Camera Tray	Polypropylene	Y	Y	Y	
125		3-3131	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
126		3-3170	One Scope 1mm-10mm	Radel	Y	Y	Y	
127		3-3180	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
128		3-3200	One Scope 1mm-10mm	Radel	Y	Y	Y	
129		9326	Super Long Scope Tray	Radel	Y	Y	Y	
130		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
131		75.0068.4	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
132	Silicone Finger Mats							
133		202025	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
134		2-1405	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
135		SP2313	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
136		SP4622	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
137		SP2719	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
138		2-3215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
139		2-1055	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
140		2-1215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
141		6-1914	Silicone Mat Bubble Bottom (no holes)	Silicone	Y	Y	Y	
142		2-4005	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
143		2111-M	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
144		2-6005	Silicone Mat-Bubble Bottom (2 piece)	Silicone	Y	Y	Y	
145		SR10	Silicone Mat-Perforated Base - Roll -Pink	Silicone	Y	Y	Y	
146								
147	Wire Baskets	75.0001.7	Wire Basket	Stainless Steel	Y	Y	Y	
148		75.005.3	Wire Basket	Stainless Steel	Y	Y	Y	
149		75.004.6	Wire Basket	Stainless Steel	Y	Y	Y	
150		75.005.9	Wire Basket	Stainless Steel	Y	Y	Y	
151		75.004.8	Wire Basket	Stainless Steel	Y	Y	Y	
152		75.004.9	Cover	Stainless Steel	Y	Y	Y	
153		75.0091.1	Cover	Stainless Steel	Y	Y	Y	
154		72.1400.5	Large Washer Basket	Stainless Steel	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	Gas-Plasma/VHP Compatibility	Additional Comments
155		72.1400.6	Large Washer Basket	Stainless Steel	Y	Y	Y	
156								
157	High-Sided Perforated Trays							
158		HDO-0952	Perforated Tray	Stainless Steel	Y	Y	Y	
159		HDO-1062	Perforated Tray	Stainless Steel	Y	Y	Y	
160		HDO-1272	Perforated Tray	Stainless Steel	Y	Y	Y	
161		HDO-1210	Perforated Tray	Stainless Steel	Y	Y	Y	
162		HDO-1610	Perforated Tray	Stainless Steel	Y	Y	Y	
163		HDO-2012	Perforated Tray	Stainless Steel	Y	Y	Y	
164								
165	Plastic Sterilization Baskets							
166		PIT2015	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
167		DINT2525	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
168		4106-H	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
169		PIT2715	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
170		MIT30251	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
171		1511-4	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
172		DT-1218	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
173		MIT40301	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
174		DT-2136	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
175		1511-7	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
176		DINT4825	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
177								
178	Perforated Stainless Steel Trays							
179		HBO-1006	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
180		20-1010	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
181		20-1015	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
182		20-1020	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
183		HBO-0515	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
184		HBO-0520	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
185		HBO-1321	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
186		HBO-1320	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
187		HBO-0624	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
188		HBO-1525	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
189		HBO-0825	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
190		HBO-0926	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
191		HBO-0630	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
192								
193	Long Scope Instrument Tray	SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	