

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
2	Laparoscopic Trays							
3		DDL-1123	Double Decker Lap Tray	Anodized Aluminum	Y	Y	Y	
4		CYSTO-2214B	Cysto Tray	Stainless Steel	Y	Y	Y	
5		DDL-4113	Poly Double Decker Tray	Radel	Y	Y	Y	
6		SDL-1130	Extra Long Lap Tray	Radel	Y	Y	Y	
7		51527	Double Decker Bariatric	Anodized Aluminum	Y	Y	Y	
8		31127 DV	Da Vinci Arm Sterilization Case	Anodized Aluminum	Y	Y	Y	
9		300605-DVC	Da Vinci Lens Tray	Anodized Aluminum	Y	Y	Y	
10		300605-DVA	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
11		300605-H	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
12		180502-INS	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
13								
14	Custom Procedure Trays							
15		30320-H	2 Level Tray System	Anodized Aluminum	Y	Y	Y	
16		31421-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
17		31123-T	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
18		31123-T1	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
19		31127-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
20								
21	Instrument Brackets for Custom Trays							
22		P-01	Post	Anodized Aluminum/Silicone	Y	Y	Y	
23		Tri-CUT BRACKET	1"universal bracket 3, 5 &10 mm	Stainless Steel/Silicone	Y	Y	Y	
24		56-UB	5" 6 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
25		710-UB	7.5" 10 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
26		1015-UB	10" 15 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
27								
28	Instru-Trays							
29		4116-H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
30		4106H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
31		41013H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
32		41116H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
33		41317H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
34		41520H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
35		HCO-1309	Perforated Tray	Stainless Steel	Y	Y	Y	
36		HCO-1510	Perforated Tray	Stainless Steel	Y	Y	Y	
37		HCO-1711	Perforated Tray	Stainless Steel	Y	Y	Y	
38		HCO-1912	Perforated Tray	Stainless Steel	Y	Y	Y	
39		HCO-1065	Perforated Tray	Stainless Steel	Y	Y	Y	
40		HCO-2116	Perforated Tray	Stainless Steel	Y	Y	Y	
41		HCO-1065-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
42		HCO-1309-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
43		HCO-1510-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
44		HCO-1711-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
45		HCO-1912-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
46								
47	Metal Micro Trays							
48		184014	Base & Cover	Anodized Aluminum	Y	Y	Y	
49		184016	Base & Cover	Anodized Aluminum	Y	Y	Y	
50		184019	Base & Cover	Anodized Aluminum	Y	Y	Y	
51		184021	Base & Cover	Anodized Aluminum	Y	Y	Y	
52		181453	Solid (no holes) Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
53		184453	Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
54		181463	Base,Cover & Holder	Anodized Aluminum	Y	Y	Y	
55		182500	Stainless Steel - 5 Clip	Stainless Steel	Y	Y	Y	
56		182062	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
57		182068	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
58								
59	Hinged Top Micro Trays							
60		30732-H	Taupe	Polypropylene	Y	Y	Y	
61		30842-H	Taupe	Polypropylene	Y	Y	Y	
62		31042-H	Taupe	Polypropylene	Y	Y	Y	
63		31055-H	White	Polypropylene	Y	Y	Y	
64		31462-H	Taupe	Polypropylene	Y	Y	Y	
65		31772-H	Blue	Polypropylene	Y	Y	Y	
66		50Z900	White	Polypropylene	Y	Y	Y	
67								
68	Hinged Top Non-Perforated Trays							
69		30732	Taupe	Propylene/Ethylene	Y	Y	Y	
70		30842	Taupe	Propylene/Ethylene	Y	Y	Y	
71		31042	Taupe	Propylene/Ethylene	Y	Y	Y	
72		31462	Taupe	Propylene/Ethylene	Y	Y	Y	
73		31772	Blue	Propylene/Ethylene				
74								
75	Specialty Micro Trays							
76		30109HN	ER Hand Trays	Anodized Aluminum	Y	Y	Y	
77		30109CC	Camera Case	Anodized Aluminum	Y	Y	Y	
78		30011KE	Key Elevators Box	Anodized Aluminum	Y	Y	Y	
79		30011CU	Currettes Box	Anodized Aluminum	Y	Y	Y	
80		30011RS	Rhoton Suction Box	Stainless Steel	Y	Y	Y	
81		30011BS	Bipolar Set Box	Stainless Steel	Y	Y	Y	
82		30011BP	Bipolar Set Box	Stainless Steel	Y	Y	Y	
83								
84	Metal Hinged Top Micro Trays							
85		31010-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
86		31010-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
87		31015-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
88		31015-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
89		31020-H	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
90								
91	Clear-Top Micro Trays							
92		0278	Base & Cover	Ultem	Y	Y	Y	
93		02781	Base & Cover	Ultem	Y	Y	Y	
94		027814	Base & Cover	Ultem	Y	Y	Y	
95		0625	Base & Cover	Ultem	Y	Y	Y	
96		062514	Base & Cover	Ultem	Y	Y	Y	
97		6540	Base & Cover	Ultem	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
98		7540	Base & Cover	Ultem	Y	Y	Y	
99		754015	Base & Cover	Ultem	Y	Y	Y	
100		100675	Single Decker Twin	Ultem	Y	Y	Y	
101		100615	Double Decker Twin	Ultem	Y	Y	Y	
102		100615B	Insert Tray for 100615	Ultem	Y	Y	Y	
103		101525	Base & Lid	Ultem	Y	Y	Y	
104		101525B	Full Insert	Ultem	Y	Y	Y	
105								
106	Silicone Plug-InSecure-It							
107		A237591K	One 4mm slot, two plugs	Silicone	Y	Y	Y	
108		A237592K	Two 4mm slots, two plugs	Silicone	Y	Y	Y	
109		A237593K	Three 4mm slots, two plugs	Silicone	Y	Y	Y	
110		A237593	Three 1mm slots, two plugs	Silicone	Y	Y	Y	
111		A237596	Six 1mm slots, two plugs	Silicone	Y	Y	Y	
112		A387596	Six 1mm slots, three plugs	Silicone	Y	Y	Y	
113		A3875912	Twelve 1mm slots, three plugs	Silicone	Y	Y	Y	
114		A387592	Two 4mm slots, three plugs	Silicone	Y	Y	Y	
115		A387593	Three 4mm slots, three plugs	Silicone	Y	Y	Y	
116		A387594	Four 4mm slots, three plugs	Silicone	Y	Y	Y	
117		A387595	Five 4mm slots, three plugs	Silicone	Y	Y	Y	
118		A2616901	One 12mm slots, two plugs	Silicone	Y	Y	Y	
119		A26216902	Two 8mm slots, two plugs	Silicone	Y	Y	Y	
120		A32759P2	One 5mm slot, three plugs	Silicone	Y	Y	Y	
121		A327593	Three 5mm slots, three plugs	Silicone	Y	Y	Y	
122		A327593L	Two 12mm slots, three plugs	Silicone	Y	Y	Y	
123		A32759P1	One 12mm slots, three plugs	Silicone	Y	Y	Y	
124								
125	Scope Trays							
126		30126	Bariatric - Single Scope Case	Anodized Aluminum	Y	Y	Y	
127		30630	Bariatric - Multi Scope Case	Anodized Aluminum	Y	Y	Y	
128		30111-SSC	Scope Tray	Stainless Steel	Y	Y	Y	
129		30118-SSC	Universal Scope Case	Stainless Steel	Y	Y	Y	
130		31042-22	Hinged Top Arthro Tray	Polypropylene	Y	Y	Y	
131		3-4110	Clear Top Arthro Tray	Radel	Y	Y	Y	
132		SP-3014	Osteotome Tray	Anodized Aluminum	Y	Y	Y	
133		3-5120	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
134		3-5130	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
135		75.0047.9M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
136		75.0083.1M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
137		75.0047.5M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
138		750083	Wire Scope Tray	Stainless Steel	Y	Y	Y	
139		31772-22	Economy Scope Tray	Polypropylene	Y	Y	Y	
140		261615	Clear Top Scope Tray	Radel	Y	Y	Y	
141		302215	Scope Tray B/L 1 set Secur-Its	Radel	Y	Y	Y	
142		31462-22	Camera Tray	Polypropylene	Y	Y	Y	
143		3-3131	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
144		3-3170	One Scope 1mm-10mm	Radel	Y	Y	Y	
145		3-3180	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
146		3-3200	One Scope 1mm-10mm	Radel	Y	Y	Y	
147		9326	Super Long Scope Tray	Radel	Y	Y	Y	
148		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
149		75.0068.4	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
150								
151	Silicone Finger Mats							

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
152		202025	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
153		2-1405	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
154		SP2313	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
155		SP4622	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
156		SP2719	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
157		2-3215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
158		2-1055	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
159		2-1215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
160		6-1914	Silicone Mat Bubble Bottom (no holes)	Silicone	Y	Y	Y	
161		2-4005	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
162		2111-M	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
163		2-6005	Silicone Mat-Bubble Bottom (2 piece)	Silicone	Y	Y	Y	
164		SR10	Silicone Mat-Perforated Base - Roll -Pink	Silicone	Y	Y	Y	
165		0278-M	Silicone Mat 2.667 x 1.5	Silicone	Y	Y	Y	
166		0625-M	Silicone Mat 2.5 W x 6L	Silicone	Y	Y	Y	
167		1006158M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
168		100615-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
169		100675-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
170		101525B-M	Finger Mat	Silicone	Y	Y	Y	
171		6540-M	Micro-Protech Mat 4 x 7.5	Silicone	Y	Y	Y	
172		7540-M	Silicone Mat	Silicone	Y	Y	Y	
173								
174	Wire Baskets							
175		75.0001.7	Wire Basket	Stainless Steel	Y	Y	Y	
176		75.005.3	Wire Basket	Stainless Steel	Y	Y	Y	
177		75.004.6	Wire Basket	Stainless Steel	Y	Y	Y	
178		75.005.9	Wire Basket	Stainless Steel	Y	Y	Y	
179		75.004.8	Wire Basket	Stainless Steel	Y	Y	Y	
180		75.004.9	Cover	Stainless Steel	Y	Y	Y	
181		75.0091.1	Cover	Stainless Steel	Y	Y	Y	
182		72.1400.5	Large Washer Basket	Stainless Steel	Y	Y	Y	
183		72.1400.6	Large Washer Basket	Stainless Steel	Y	Y	Y	
184								
185	High-Sided Perforated Trays							
186		HDO-0952	Perforated Tray	Stainless Steel	Y	Y	Y	
187		HDO-1062	Perforated Tray	Stainless Steel	Y	Y	Y	
188		HDO-1272	Perforated Tray	Stainless Steel	Y	Y	Y	
189		HDO-1210	Perforated Tray	Stainless Steel	Y	Y	Y	
190		HDO-1610	Perforated Tray	Stainless Steel	Y	Y	Y	
191		HDO-2012	Perforated Tray	Stainless Steel	Y	Y	Y	
192								
193	High-Sided Non-Perforated Trays							
194		HDO-1062-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
195		HDO-1210-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
196		HDO-1272-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
197		HDO-1610-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
198		HDO-2012-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
199		HDO-952-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
200								
201	Plastic Sterilization Baskets							
202		PIT2015	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
203		DINT2525	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
204		4106-H	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
205		PIT2715	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
206		MIT30251	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
207		1511-4	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
208		DT-1218	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
209		MIT40301	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
210		DT-2136	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
211		1511-7	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
212		DINT4825	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
213								
214	Perforated Stainless Steel Trays							
215		PSST-009	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
216		PSST-007	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
217		PSST-005	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
218		PSST-010	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
219		PSST-001	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
220		PSST-002	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
221		PSST-012	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
222		PSST-003	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
223		PPST-011	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
224		PSST-004	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
225		PSST-013	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
226		PSST-008	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
227		PSST-006	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
228								
229	Long Scope Instrument Tray							
230		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
231								
232	Sterilizing Basket							
233		PIT3025	Sterilizing Basket 11.82 x 9.85 x 2.05H	Stainless Steel	Y	Y	Y	
234		PIT5739	Sterilizing Basket 21.81 x 15.44 x 2.05H	Stainless Steel	Y	Y	Y	
235								
236	Nic-Nac Washer Baskets							
237		CIL-004	Fine Mesh Basket with Snap Closure	Stainless Steel	Y	Y	Y	