



Instructions for Use: Multi-purpose Labels

Brand Name of Product	Multi-purpose Labels
Generic Name of Product	Labels
Product Code Number(s)	HM-CID-613, HM-52483-HLD, CLN-4X4
Purpose of Product	Provides the flexibility to communicate clean or biohazard supplies on a label designed specifically for your department.
Range of Applications for Product	<ul style="list-style-type: none"> To be placed on equipment, trays or anywhere key communication is required.
Key specifications of product	<ul style="list-style-type: none"> “Clean” Labels “In Use” “Dirty” Labels Adhesive back

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	Store in room temperature.
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	Peel off adhesive back and place firmly on equipment, instrument trays or desired location.
Preparation for Use	Peel adhesive from the back. Apply to clean, dry, room temperature surfaces.
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Peel adhesive from the back of the label 2. Place on desired location 3. Use two fingers to press firmly from the middle and move your fingers outward to the edges of the label
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	Do not use on items that are going to be sterilized. Labels are not autoclavable.
Disposal	Single use, treat as biohazard material. Dispose of according to your facilities guidelines.

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	
Storage:	

Additional Information:	Do not use on items that are going to be sterilized.
Related Healthmark Products	
Other Product Support Documents	
Reference Documents	Label Price List, Label Brochure
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-09-29 Ralph J Basile