

Instructions for Use: Clean & Dirty Reversible Signs

Brand Name of Product	Clean & Dirty Reversible Signs
Generic Name of Product	Label Holders
Product Code Number(s)	2024-4060, 2022-7555, 2023-1058, CD-2022, CD-2023, CD-2024
Intended Use	A secure method for displaying “Clean” or “Dirty” on nearly any cart or trolley to communicate the status of transported instruments and equipment.
Range of Applications for Product	Protects signs, photographs, business cards, labels, diagrams and other documents.
Key specifications of product	<ul style="list-style-type: none"> • Magnet loops • Vinyl sheet • Reversible sign displaying “Clean” & “Dirty”

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	When the signs are not in use, please DO NOT stack more than 10 signs together. This will protect the quality of the product between uses.
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place loops around the handle of the cart, trolley or device. 2. Secure loops with the connecting magnets. 3. Use “Clean” side when transporting clean instruments and equipment. 4. Flip the sign to “Dirty” when transporting used or dirty instruments or equipment.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	When the signs are not in use, please DO NOT stack more than 10 signs together. This will protect the quality of the product between uses.
Disposal	

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	Clean with warm water and soap.
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	
Storage:	
Additional Information:	Do NOT put in the cart washer, automated washer or sterilizer. It will distort the vinyl.

Related Healthmark Products	
Other Product Support Documents	Bins Galore Price List, Bins Galore Brochure
Reference Documents	

Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com
----------------------------------	---

2018-01-16 Matt Smith