


Instructions for Use: Clip-On Basket Tag

Brand Name of Product	Clip-On Basket Tag
Generic Name of Product	Identification Tags
Product Code Number(s)	COBT OR, CBTOR-1, CBTOR-25, COBT RD, CBTRD-1, CBTRD-25, COBT BL, CBTBL-1, CBTBL-25, COBT YL, CBTYL-1, CBTYL-25, COBT GN, CBTGN-1, CBTGN-25, COBT PK, CBTPK-1, COBT BK, CBTBK-1, CBTBK-25, CBTPK-25, COBT PP, CBTPP-1, CBTPP-25
Purpose of Product	These labels facilitate specific markings and remarks for trays. They are designed to minimize misunderstandings about the contents of a tray and to notify about needed maintenance of contents.
Range of Applications for Product	Mesh baskets
Key Specifications of Product	<ul style="list-style-type: none"> • Clips on the back ideal for attaching to mesh baskets • Heat resistant up to 135 °C or 275 °F • Available in eight different colors

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	Sold in packs of 25
Shelf Life	

Instructions for Using Product	
Description of Use(s)	The basket tag is heat, water and chemical resistant and thus can stay with the basket all the way through the process.
Preparation for Use	
Diagrams (drawings, pictures)	
Steps for Use of Product	Attach the clips of the tag to the basket.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	Do not handle basket tags until they cool down to room temperature.
Disposal	

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	Compatible with steam, EO and VH ₂ O ₂ sterilization.
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	Label Price List, Label Brochure
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2021-03-29 Suzanne Latta