



Brand Name of Product	Custom Labels
Generic Name of Product	Custom Labels
Product Code Number(s)	CSW-116, CSW-254, CSW-289, CSW-302, CSW-GS, CSW-IMP, CSW-RCS, CSW-RRS, CSW-SPU, CSW-SSS, DC-100214-71-BLUE, CAL-33, CAL-8722, CAL-873
Intended Use	Custom labels with specific messages designed to fit individual department needs
Range of Applications for Product	Wide range depending on message intended to communicate
Key specifications of product	<ul style="list-style-type: none"> • Various sizes • Various colors • Steam, Gas-Plasma and EO autoclavable optional • CAL labels available up to 1000. • CSW and DC-100214-71-BLUE require min order (see price list)

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	Custom labels with specific messages designed to fit individual department needs
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Order labels by completing order sheet. Select whether adhesive needs to be autoclavable. (CAL labels are always steam sterilizable) 2. Peel adhesive from the back of the label. 3. Place on desired location. 4. Use two fingers to press firmly from the middle and move your fingers outward to the edges of the label.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	<ul style="list-style-type: none"> • Steam autoclavable labels will also withstand cleaning in a mechanical washer.
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	<ul style="list-style-type: none"> • If subjected to sterilization, be sure to select a label compatible with the sterilization method to be used. • CAL labels are compatible with steam sterilization cycles: (270°F – for 4 minutes / 275°F for 3 minutes.
Storage:	
Additional Information:	

Related Healthmark Products	
Other Product Support Documents	Labeling Brochure, Labeling Price List

Reference Documents	
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-09-15 msmith