



**Instructions for Use: Custom Labels**

<b>Brand Name of Product</b>	Custom Labels
<b>Generic Name of Product</b>	Custom Labels
<b>Product Code Number(s)</b>	CSW-116, CSW-254, CSW-289, CSW-302, CSW-GS, CSW-IMP, CSW-RCS, CSW-RRS, CSW-SPU, CSW-SSS, DC-100214-71-BLUE, CAL-33, CAL-8722, CAL-873
<b>Intended Use</b>	Custom labels with specific messages designed to fit individual department needs
<b>Range of Applications for Product</b>	Wide range depending on message intended to communicate
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Various sizes</li> <li>• Various colors</li> <li>• Steam, Gas-Plasma and EO autoclavable optional</li> <li>• CAL labels available up to 1000.</li> <li>• CSW and DC-100214-71-BLUE require min order (see price list)</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	Custom labels with specific messages designed to fit individual department needs
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Order labels by completing order sheet. Select whether adhesive needs to be autoclavable. (CAL labels are always steam sterilizable)</li> <li>2. Peel adhesive from the back of the label.</li> <li>3. Place on desired location.</li> <li>4. Use two fingers to press firmly from the middle and move your fingers outward to the edges of the label.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	<ul style="list-style-type: none"> <li>• Steam autoclavable labels will also withstand cleaning in a mechanical washer.</li> </ul>
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• If subjected to sterilization, be sure to select a label compatible with the sterilization method to be used.</li> <li>• CAL labels are compatible with steam sterilization cycles: (270°F for 4 minutes / 275°F for 3 minutes.</li> </ul>
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Labeling Brochure, Labeling Price List

<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>