




**Instructions for Use: Custom Vinyl Labels**

<b>Brand Name of Product</b>	Custom Vinyl Labels
<b>Generic Name of Product</b>	Labels
<b>Product Code Number(s)</b>	CVL14, CVL16, CVL18, CVL112, CVL24, CVL26, CVL28, CVL212, CVL34, CVL36, CVL38, CVL312, CVL44, CVL46, CVL48, CVL412, CVL55, CVL56, CVL58, CVL512, CVL66, CVL68, CVL612, CVL76, CVL77, CVL78, CVL712, CVL86, CVL88, CVL812, CVL96, CVL98, CVL99, CVL912
<b>Purpose of Product</b>	Provides the flexibility to communicate any message on a label designed specifically for your department.
<b>Range of Applications for Product</b>	<ul style="list-style-type: none"> <li>To be placed on equipment, trays or anywhere key communication is required.</li> </ul>
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Material is vinyl</li> <li>Sizes range from 1”x 4” to 4” x 12”</li> <li>Custom sizes available for quotation</li> <li>Comes in a variety of colors</li> <li>Adhesive back</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Peel off adhesive back and place firmly on equipment, instrument trays or desired location.		
<b>Preparation for Use</b>	Peel off adhesive from the back.		
<b>Diagrams (drawings, pictures)</b>			
	1	2	3
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Peel adhesive from the back of the label.</li> <li>2. Place on desired location.</li> <li>3. Use two fingers to press firmly from the middle and move your fingers outward to the edges of the label.</li> </ol>		
<b>Interpretation of Results</b>			
<b>Contraindications of Test Results</b>			
<b>Documentation</b>			
<b>Special Warnings and Cautions</b>	Do not use on items that are going to be sterilized. Labels are not autoclavable.		
<b>Disposal</b>			

**Reprocessing Instructions**

<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	Wash with any mild detergent. Do not exceed temperature of 180°F.
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	

<b>Storage</b>	
<b>Additional Information</b>	Do not use on items that are going to be sterilized.
<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	
<b>Reference Documents</b>	Label Price List, Label Brochure
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2021-02-08 Ralph J Basile