




Instructions for Use: Custom Magnetic Labels

Brand Name of Product	Magnetic Labels
Generic Name of Product	Labels
Product Code Number(s)	CVL-M14, CVL-M16, CVL-M18, CVL-M112, CVL-M24, CVL-M26, CVL-M28, CVL-M212, CVL-M34, CVL-M36, CVL-M38, CVL-M312, CVL-M44, CVL-M46, CVL-M48, CVL-M412, CVL-M55, CVL-M56, CVL-M58, CVL-M512, CVL-M76, CVL-M78, CVL-M712, CVL-M86, CVL-M88, CVL-M812, CVL-M96, CVL-M98, CVL-M99, CVL-M912
Purpose of Product	Magnetic labels provide solution for removing and re-applying labels to metal surfaces. Use them to communicate important messages on cabinets, sinks or any other metal surface in your department.
Range of Applications for Product	<ul style="list-style-type: none"> To be placed on metal equipment, trays or anywhere key communication is required.
Key specifications of product	<ul style="list-style-type: none"> Material is vinyl Sizes range from 1”x 4” to 4” x 12” Custom sizes available for quotation Comes in a variety of colors Magnetic back

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product

Description of Use (s)	Magnetic back. Place firmly on equipment, instrument trays or desired location.
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> Find desired metal surface. Stick magnetic side of label onto surface. Make sure label is secure.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	Do not use on items that are going to be sterilized. Labels are not autoclavable.
Disposal	

Reprocessing Instructions

Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	Wash with any mild detergent. Do not exceed temperature of 180 ° F
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	
Storage:	
Additional Information:	Do not use on items that are going to be sterilized.

Related Healthmark Products	
Other Product Support Documents	
Reference Documents	Label Price List, Label Brochure
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-12-06 Ralph J Basile