


## Instructions for Use: Stainless Steel Tray Tags

<b>Brand Name of Product</b>	Stainless Steel Tray Tags
<b>Generic Name of Product</b>	Tray Tags
<b>Product Code Number(s)</b>	01444, 01452, 01454
<b>Purpose of Product</b>	Alternative for identifying trays, containers, and other key assets.
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Stainless steel</li> <li>Three styles to choose from</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Con</b>	N/A
<b>Shelf Life</b>	N/A

Instructions for Using Product	
<b>Description of Use(s)</b>	Attach tray tags to trays for identification.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	 <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <span>Clip-On Style</span> <span>Round Style</span> <span>Tie-In Rectangular Style</span> </div>
<b>Steps for Use of Product</b>	Three styles for use: <ol style="list-style-type: none"> <li>Clip-on will snap on the rim of the tray.</li> <li>Tie-in round and rectangle tags have pre-cut holes to use ties tags for mounting to trays and other equipment.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

Reprocessing Instructions	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	May be processed with the basket of instruments in accordance with the cleaning procedure for those instruments.
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>Material compatible with steam sterilization at 132°C (270°F) for 4 minutes, 135°C (275°F) for 3 minutes.</li> <li>Compatible with EO and VH2O2 Sterilization.</li> </ul>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	Labeling Brochure, Labeling Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2021-04-22 Suzanne Latta