



04/05/2010

CLINIPAK TRAY TAGS – COMPATIBLE STERILISATION METHODS

I hereby confirm that Clinipak Tray Tags with the Following Product Codes are compatible with Steam (up to 137°C), Ethylene Oxide and Plasma Sterilisation Methods:

TAG001 – TAG026 inclusive
TAG55WH

A handwritten signature in black ink that reads "Andy Wigglesworth". The signature is written in a cursive style with a horizontal line underneath.

Andy Wigglesworth
QA Manager
Clinipak Ltd

Tel: 01628 648091
Email: qualitycontrol@clinipak.co.uk
Fax 01628 810760

Clinipak Limited · Beech House · Knaves Beech Business Centre · Loudwater · High Wycombe · Bucks · HP10 9SD
Tel: +44 (0) 1628 810626 Fax: +44 (0) 1628 810760 Email: sales@clinipak.co.uk www.clinipak.co.uk