






**Instructions for Use: Custom Vinyl Labels**

<b>Brand Name of Product</b>	Custom Vinyl Labels
<b>Generic Name of Product</b>	Labels
<b>Product Code Number(s)</b>	CVL14, CVL16, CVL18, CVL112, CVL24, CVL26, CVL28, CVL212, CVL34, CVL34, CVL36, CVL38, CVL312, CVL44, CVL46, CVL48, CVL412
<b>Purpose of Product</b>	Provides the flexibility to communicate any message on a label designed specifically for your department.
<b>Range of Applications for Product</b>	<ul style="list-style-type: none"> <li>To be placed on equipment, trays or anywhere key communication is required.</li> </ul>
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>Material is vinyl</li> <li>Sizes range from 1”x 4” to 4” x 12”</li> <li>Custom sizes available for quotation</li> <li>Comes in a variety of colors</li> <li>Adhesive back</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Peel of adhesive back and place firmly on equipment, instrument trays or desired location.		
<b>Preparation for Use</b>	Peel adhesive from the back		
<b>Diagrams (drawings, pictures):</b>			
	1	2	3
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Peel adhesive from the back of the label</li> <li>2. Place on desired location</li> <li>3. Use two fingers to press firmly from the middle and move your fingers outward to the edges of the label</li> </ol>		
<b>Interpretation of Results</b>			
<b>Contraindications of Test Results</b>			
<b>Documentation</b>			
<b>Special Warnings and Cautions</b>	Do not use on items that are going to be sterilized. Labels are not autoclavable.		
<b>Disposal</b>			

**Reprocessing Instructions**

<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	Wash in accordance of product IFU that the label is being placed on.
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	Do not use on items that are going to be sterilized.

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	
<b>Reference Documents</b>	Label Price List, Label Brochure
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com

2013-08-27 Ralph J Basile