




Healthmark Industries Company

LABELING & IDENTITY

A complete line of inventory control, message and customs labels for identification, tracking and communication.




Instrument Marking
Page 2



Tray Marking & Tracking
Page 3




Label Gun & Labels
Page 4



Color Code Label Tape & Dots
Page 5



Custom Autoclave Tape
Page 6



Message Labels
Page 6

SURGERY LINE-UP				4-19-08
O.R.	TIME	SURGEON	STATUS	PROCEDURE
	7:30A	MARKS		
	9:00A	Jones		
	7:30A	STILLS		
	9:30A	HANSON		
	7:00A	KENT		
	8:15A	CARLSON		
	10:00A	KENT		
	8:00A	MASON		
	10:00A	MARKS		

OR. 1	OR. 3
OR. 2	OR. 4

Instant Grids
Page 9



Signs
Page 10



Item Tracing
Page 12



INSTRUMENT & TRAY TRACKING



Tray Identification Tags


Resin bonded tags designed for tray identification. Tags may be used with bar code labels or hand written notes. All instruments remain with the sterilization basket throughout the collection, decontamination and sterilization process. X-ray detectable.

Code No	Description	Width	Length
RIN001	Stainless Ring, 100 ea.	1.00	
TAG001	White, 100 ea.	3.50	2.20
TAG011	White, 100 ea.	3.50	1.50
TAG012	Tan, 100 ea.	3.50	1.50
TAG013	Yellow, 100 ea.	3.50	1.50
TAG014	Orange, 100 ea.	3.50	1.50
TAG015	Red, 100 ea.	3.50	1.50
TAG016	Blue, 100 ea.	3.50	1.50
TAG017	Green, 100 ea.	3.50	1.50
TAG018	Grey, 100 ea.	3.50	1.50
TAG019	Purple, 100 ea.	3.50	1.50
TAG020	Black, 100 ea.	3.50	1.50
TAG021	Lime, 100 ea.	3.50	1.50
TAG022	Pink, 100 ea.	3.50	1.50
TAG023	Turquoise, 100 ea.	3.50	1.50
TAG024	Mustard, 100 ea.	3.50	1.50
TAG025	Terracotta, 100 ea.	3.50	1.50
TAG026	Royal Blue, 100 ea.	3.50	1.50
TIE001	Plastic Tie, 100 ea.		
TIE002	Plastic Tie, 100 ea.		

Convenient Instrument Marking Tape

Instrument marking tape comes in a convenient binder which is indexed by color-coded separators. Tape is precut in two lengths, 1^{1/8}"L x 3/8"W (116 per sheet) and 2^{1/4}"L x 3/8"W (18 per sheet) for easy application. Steam autoclavable. Available in many colors. Check the Labels Products Price List for a complete listing.

Code No	Description
MKT003	Binder w/ index & 1 sheet ea. Zebra
MKT001	Binder w/ index
MKT002	Binder w/ index & 1 sheet. Plain
MKT004	Complete Binder Set
(see below)	Refill Packs of 5 Sheets Plain
(see below)	Refill Packs of 5 Sheets Zebra




INSTRUMENT & TRAY TRACKING



Code No	Description	Length	Width
403005	White Tray Label	6.50	4.00
403115	White Tray Label	9.50	2.00
403125	Preprinted Tray Label	9.50	2.00
403225	White Tray Label	8.50	2.50
403225 BLUE	Blue Tray Label	8.50	2.50
403225 GREEN	Green Tray Label	8.50	2.50
403225 HANG TIME	Label with Month and Dates	8.50	2.50
403225 LOANER	"LOANER" White Tray Label	8.50	2.50
403225 PINK	Pink Tray Label	8.50	2.50
403225 RED	Red Tray Label	8.50	2.50
403225 SLOTTED	White, Slotted for Count Sheet	8.50	2.50
403225 SLOTTED PINK	Pink, Slotted for Count Sheet	8.50	2.50
403225 YELLOW	Yellow Tray Label	8.50	2.50

Improve Tracking with the Instrument Tray Labels

These labels facilitate specific markings and remarks for sterile or contaminated instrument trays. The PP-based label is suitable for sterilization and are resistant to washing and disinfection and thus can stay with a tray of instruments all the way through the decontamination process. Packaged 500 per box. Need a special design? Call for a custom label quote.



Stainless Steel Tray Tags

Available in 3 styles, these stainless steel tags are a great alternative for identifying trays, containers and other key assets. Choose from clip-on, tie-on in round and rectangle.

Code No	Description	Length	Width
01444	Stainless Steel Tag	3.50	1.50
01452	Stainless Steel Badge	0.75	0.00
01454	Clip-on Stainless Steel	3.50	1.50



LABEL GUN SYSTEM

Hand-held Labeling Systems



**Free Gun!
Ask for
Details**

Healthmark's Hand-held Labeling System has been designed specifically with Central Sterile in mind for the labeling of any items that you need to be autoclaved - from sterile linen packs to instrument trays, and peel pouches. It also has great versatility for labeling in other hospital departments such as dietary and pharmacy.

This system gives you the flexibility to print two lines of information with up to eight characters per line. The available characters include two gregorian dates for printing expiration and autoclave dating with additional bands for load and sterilizer numbers.

NEW! Gas-Plasma compatible labels are available. "S" at end of code number is the Gas-Plasma Compatible label.

Code No.	Description	Pkg
H-362026	HANDHELD LABEL GUN	EACH
H369020	EVENT RELATED LABEL - RED (steam & EtO)	CS-12
H369021	EVENT RELATED LABEL - BLUE (steam & EtO)	CS-12
H369023	EVENT RELATED LABEL - BLACK (steam & EtO)	CS-12
H369024	EVENT RELATED LABEL - GREEN (steam & EtO)	CS-12
H369010	EXPIRATION LABEL - BLUE (steam & EtO)	CS-12
H369012	EXPIRATION LABEL - RED (steam & EtO)	CS-12
H369013	EXPIRATION LABEL - BLACK (steam & EtO)	CS-12
H369014	EXPIRATION LABEL - GREEN (steam & EtO)	CS-12
H369020S	EVENT RELATED STERRAD LABELS - RED (Sterrad, steam & EtO)	CS-12
H369021S	EVENT RELATED STERRAD LABELS - BLUE (Sterrad, steam & EtO)	CS-12
H369023S	EVENT RELATED STERRAD LABELS - BLACK (Sterrad, steam & EtO)	CS-12
H369024S	EVENT RELATED STERRAD LABELS - GREEN (Sterrad, steam & EtO)	CS-12
H369010S	EXPIRATION LABELS STERRAD - BLUE (Sterrad, steam & EtO)	CS-12
H369012S	EXPIRATION LABELS STERRAD - RED (Sterrad, steam & EtO)	CS-12
H369013S	EXPIRATION LABELS STERRAD - BLACK (Sterrad, steam & EtO)	CS-12
H369014S	EXPIRATION LABELS STERRAD - GREEN (Sterrad, steam & EtO)	CS-12

Event Related Labels

Green




H369024/H369024S

Red



H369020/H369020S

Black



H369023/H369023S

Blue



H369021/H369021S


Expiration Labels

Green



H369014/H369014S
(steam & eto / gas plasma)

Red



H369012/H369012S

Blue




H369010/H369010S

Black



H369013/H369013S

STERILIZATION & IDENTITY



Color Code Label Tape


This tape readily accepts ballpoint, pencil or permanent markers. With a service temperature of -70°F to +250°F, it can be used in freezer, refrigerator, autoclave or incubator. This tape will stay put until you purposely remove it and comes off cleanly with no adhesive transfer. Available in 2 convenient length rolls, 5 widths and 26 colors.

Code No.	Length	Width	Pkg
LT-512	500	0.50	ROLL
LT-534	500	0.75	ROLL
LT-510	500	1.00	ROLL
LT-614	2160	0.25	ROLL
LT-612	2160	0.50	ROLL
LT-610	2160	1.00	ROLL
LT-620	2160	2.00	ROLL



Color Dots

Dating, coding, categorizing is made simple with these Color Dots. Supplied in sheet form or on rolls form in 7 sizes and 26 colors.



Dia-meter	Sheet Code No.	Roll Code No.	Pkg
0.25	LDL-025	LDR-025	PKG-1000
0.50	LDL-050	LDR-050	PKG-1000
0.75	LDL-075	LDR-075	PKG-1000
1.00	LDL-100	LDR-100	PKG-1000
1.25	LDL-125	LDR-125	PKG-1000
1.50	LDL-150	LDR-150	PKG-1000
2.00	LDL-200	LDR-200	PKG-1000



Removable Page Marker Tags

- Repositionable with restickable adhesive.
- Half-adhesive leaves half of the tag exposed.
- Color code for different people or projects.

Code No.	Description
RT01	REDI-TAG - YELLOW
RT02	REDI-TAG - RED
RT03	REDI-TAG - BLUE
RT04	REDI-TAG - DARK BLUE
RT05	REDI-TAG - GOLD
RT06	REDI-TAG - SILVER
RT07	REDI-TAG - ORANGE
RT08	REDI-TAG - GREEN
RT09	REDI-TAG - DARK GREEN
RT10	REDI-TAG - PINK
RT12	REDI-TAG - BROWN
RT13	REDI-TAG - BLACK
RT14	REDI-TAG - PURPLE
RT18	REDI-TAG - MINT
RT21	REDI-TAG - FLUORESCENT RED
RT22	REDI-TAG - FLUORESCENT PINK



AUTOCLAVE TAPE & MESSAGE LABELS

Custom Imprinted Autoclave Tape

Imprinted tape provides easy content identification. Over 200 standard imprints available. Specify your own custom imprint at the same low price (8 roll minimum) or select from 6 Color Code tapes.



Code No.	Description
ATP-115	AUTOCLAVE TAPE "BASIN SETS"
ATP-115-04	AUTOCLAVE TAPE "BASIN SETS" PRINTED IN RED
ATP-130	AUTOCLAVE TAPE "CIRCUMCISION TRAY"
ATP-130-02	CUSTOM AUTOCLAVABLE TAPE-YELLOW "CIRCUMCISION"
ATP-130-03	AUTOCLAVE TAPE "CIRCUMCISION TRAY" GREEN INK
ATP-130-4	CUSTOM AUTOCLAVABLE TAPE RED "CIRCUMCISION"
ATP-154-4	"DRESSING TRAY" AUTOCLAVABLE TAPE -RED
ATP-223	AUTOCLAVE TAPE "O.B. INSTRUMENTS"
ATP-223-4	CUSTOM AUTO. TAPE W/ RED INK "O.B. INSTRUMENTS"
ATP-239	AUTOCLAVE TAPE "PREP BASINS"
ATP-250-4	"SCISSORS" AUTOCLABLE TAPE - RED
ATP-261	AUTOCLAVE TAPE "SUTURE REMOVAL SET"
ATP-262	AUTOCLAVE TAPE "SUTURE SETS"
ATP-269	AUTOCLAVE TAPE "THORACTOMY"
ATP-269-23	AUTOCLAVE TAPE "THORACOTOMY TRAY" BROWN INK
ATP-276	AUTOCLAVE TAPE "TOWELS"
ATP-634	PRINTED AUTOCLAVE TAPE 3/4" X 60 YDS
ATP-634-4	CUSTOM PRINTED AUTOCLAVE TAPE - RED
ATP-634C-04	AUTOCLAVE TAPE RED
ATP-634C-23	AUTOCLAVE TAPE BROWN
ATP-634C-24	AUTOCLAVE TAPE DARK BLUE
ATP-634C-26	CUSTOM AUTOCLAVABLE TAPE

PREPRINTED MESSAGE LABELS


actual size



CSW-xxx



CSW-GW




CSW-113



CSW-116



CSW-SSS



CSW-165

Call:800-521-6224 • Fax:586-491-2113 • www.HMARK.com

PREPRINTED MESSAGE LABELS



PREPRINTED MESSAGE LABELS

actual size

**WARNING
GAS STERILIZE
ONLY**

CSW-1212

STERILIZED
Indefinite Shelf Life
Unless Damaged or Opened

CSW-1216

AUTOCLAVED ON _____
EXPIRES ON _____

RETURN TO CENTRAL SERVICE
FOR RESTERILIZATION

CSW-1215

**This Tray is
INCOMPLETE**

Missing _____

_____ Date _____

CSW-302

**STERILITY IS GUARANTEED
UNTIL THE PACKAGE IS
DAMAGED OR OPENED**

CSW-254

STERILIZED
JANUARY _____ LOAD _____

STERILIZED
FEBRUARY _____ LOAD _____

STERILIZED
MARCH _____ LOAD _____

STERILIZED
APRIL _____ LOAD _____

STERILIZED
MAY _____ LOAD _____

STERILIZED
JUNE _____ LOAD _____

STERILIZED
JULY _____ LOAD _____

STERILIZED
AUGUST _____ LOAD _____

STERILIZED
SEPTEMBER _____ LOAD _____

STERILIZED
OCTOBER _____ LOAD _____

CSW-xxx
Specify Month

Code No.	Description	Length	Width	Pkg
CSW-116	EXPIRATION DATE	2.00	0.88	ROLL
CSW-1216	"STERILIZED INDEFINITE SHELF LIFE"	2.50	1.00	RL-1000
CSW-254	STERILITY IS GUARANTEED	2.25	0.88	RL-500
CSW-289	"PLEASE IMMERSE USED INSTRUMENTS IN DETERGENT..."	3.00	0.88	ROLL
CSW-302	THIS TRAY IS INCOMPLETE LABEL - ORANGE	2.25	1.75	ROLL
CSW-GS	GAS STERILIZED	2.50	0.88	ROLL
CSW-IMP	"IMPLANTABLE" LABEL BLACK ON RED	3.00	0.88	RL-320
CSW-RCS	"RETURN TO CENTRAL SUPPLY" LABEL	2.25	0.88	ROLL
CSW-RRS	REUSABLE RETURN TO SPD	2.25	0.88	ROLL
CSW-SPU	"SINGLE PATIENT USE ONLY-THEN..."	2.25	0.88	ROLL
CSW-SSS	STERRAD SYSTEM STERILIZED	2.25	0.88	ROLL

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com



CUSTOM LABELS



Design Your Own Labels





With a few easy steps, you can design a label of your own.

- Step 1 - Pick the Adhesive.
- Step 2 - Choose from these standard Label Sizes (call for custom sizes).
- Step 3 - Select from 24 Material Colors.
- Step 4 - Select from 24 Ink Colors.
- Step 5 - Sketch your layout within the appropriate label size (below).
- Step 6 -Fax this form to 586-774-6473. Indicate the quantity desired and voila, you have your own label!

1 ADHESIVES

- Permanent
- Removable
- Ultra Removable
- STERRAD Compatible

2 SIZE



3

MATERIAL COLOR


4

TEXT COLOR

#2-Yellow	#15-Lime
#3-Green	#18-Aqua
#6-Red	#17-Chartruse
#5-Orange	#18-Gold
#8-Blue	#19-Silver
#7-Pink	#20-Aspen
#9-Tan	#21-Dark Red
#8-Salmon	#22-Maroon
#10-Copper	#23-Dark Brown
#11-Rose	#24-Dark Blue
#12-Violet	#25-Dark Green
#13-Lavender	#26-Black
#14-Gray	#1-White

5 LAYOUT

INSTA GRIDS



INSTA GRIDS

Insta Grids are grids that are ready to use in seconds! They are available for use on a magnetic board or peel & stick version for use on any whiteboard. Just stick on the board creating an instant grid. Each square is perfectly shaped and presents a professional appearance. Insta Grids can be customized to any application including creating larger rectangles by cutting grid cells to size. Connect grids together to add columns and rows as needed. Insta Grids are available in six colors. Consider ordering different colors to color code your data. One kit contains six 15" x 23" grid panels to complete a 4' x 4' grid board.

Magnetic Style

Code No.	Color	Width	Height
850	Red	15.33	23.00
851	Orange	15.33	23.00
852	Yellow	15.33	23.00
853	Green	15.33	23.00
854	Blue	15.33	23.00
855	Black	15.33	23.00

MAGNA TAGS

Quickly label any magnetic surface with the Magna Tags. Attach your own label on the surface of the tag, or use a choose from six different colors and write directly on the erasable surface with a magic marker.



Peel & Stick Style

Code No.	Color	Width	Height
860	Red	15.33	23.00
861	Orange	15.33	23.00
862	Yellow	15.33	23.00
863	Green	15.33	23.00
864	Blue	15.33	23.00
865	Black	15.33	23.00

Magnet Tags

Code No.	Color	Width	Height
870	Red	3.50	1.25
871	Orange	3.50	1.25
872	Yellow	3.50	1.25
873	Green	3.50	1.25
874	Blue	3.50	1.25
875	Black	3.50	1.25



FINE SIGNS

FINE SIGNS

Select from our stock signs or design your own (see next page) to alert staff and/or visitors to important information.



FS1624-TO



FS2412-WC



FS1624-MSK



FS2424-FLAM



FS1224-AIDE



FS2406-CATH



FS1624-BIO



FS2412-EYE



FS2406-P&P



FS1224-LIFT



FS2424-KAC



FS2424-STOP



FS2412-CEILING



FS2412-PPE



FS2424-RSTR



FS2424-LABCOAT



Design Your Own Sign

With a few easy steps, you can design a sign of your own.

- Step 1 - Choose from these standard sign sizes (call for custom sizes).
- Step 2 - Select from 6 Background (poster) Colors.
- Step 3 - Select from 6 Foreground (vinyl) Colors.
- Step 5 - Sketch your layout within the appropriate guides for you size poster (below).
- Step 6 - Fax this form to 586-774-6473. Indicate the quantity desired and voila, you have your own sign!

1 SIZE

Select Width:

- 6"
- 12"
- 16"
- 24"

Select Height:

- 6"
- 12"
- 16"
- 24"

FS0606	6" wide x 6" high
FS0612	6" wide x 12" high
FS0616	6" wide x 16" high
FS0624	6" wide x 24" high
FS1206	12" wide x 6" high
FS1212	12" wide x 12" high
FS1216	12" wide x 16" high
FS1224	12" wide x 24" high

FS1606	16" wide x 6" high
FS1612	16" wide x 12" high
FS1616	16" wide x 16" high
FS1624	16" wide x 24" high
FS2406	24" wide x 6" high
FS2412	24" wide x 12" high
FS2416	24" wide x 16" high
FS2424	24" wide x 24" high


3 BACKGROUND COLOR



4 FOREGROUND COLOR



5 LAYOUT





TRACE MARKS

TRACE MARKS

Create a place for everything and put everything in its place. Trace



Marks are the quick, easy and professional looking method for identifying where instruments and other items should be replaced after use. Missing



items can also be quickly recognized and found. Use in carts, on peg boards, in maintenance areas, work stations or as part of a quality inspection process. Trace Marks are incredibly easy to use. Just trace out the instrument, cut-out your pattern and the adhere in place. Available in six different colors, which you can use to color code applications.



Code No.	Color	Width	Length
403	Black	12.25	180.00
404	White	12.25	180.00
405	Green	12.25	180.00
406	Red	12.25	180.00
407	Orange	12.25	180.00
408	Blue	12.25	180.00

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com