

Date: 6/26/2018

Healthmark Industries introduces the Foam Tipped Swab to its ProSys™
Instrument Care product line

Foam Tipped Swab



Healthmark Industries is delighted to introduce the Foam Tipped Swab to its ProSys™ Instrument Care product line. Designed for cleaning equipment, the swab consists of a polyurethane foam tip with a polypropylene handle. The foam tipped applicator is produced with reticulated 100 ppi medical grade foam, is 13.49 mm in diameter and has a length of 25.4 mm. It is thermally bonded to a rigid rectangular polypropylene shaft that has a diameter of 3.18 mm, width of 4.95 mm and length of 127mm. The swab is non-linting and residue free, as well as soft and non-abrasive. It can be used with both IPA and DI water, or dry and can be purchased in a quantity of 50.

Part#	Description	Pkg	Price
FTS-001	5 Inch Foam Tipped Swab with Polypropylene Handle	BX-50	\$100.00

More information available upon request.