

Date: 4/26/2019

Healthmark Industries introduces the Glove Box Dispenser to its Personal Protection Equipment product line

Glove Box Dispenser for Xtra-Long



Healthmark Industries is delighted to introduce the Glove Box Dispenser for Xtra Long Gloves (16”) to its Personal Protection Equipment product line. Manufactured from clear PETG Plastic-Copolyester resin, the Glove Box Dispenser is designed for boxed disposable glove storage, made to withstand heavy-duty applications and handling, as well as reduce losses due to breakage. The Glove Box Dispensers are offered in two different styles, a 0.8 lbs. single dispenser measuring 5.45 W x 16.10 H x 2.95 D Inches that accommodates one box of extended cuff gloves, and a 1.8 lbs. triple dispenser with dividers measuring 15.70 W x 16.10 H x 2.70 D Inches that holds three boxes of extended cuff gloves. The clear design matches any décor and enables an unobstructed view of glove package information. The Glove Box Dispenser is equipped with keyholes for wall mounting, which saves counter space while keeping gloves within easy reach. Both styles are available for individual purchase.

Part#	Description	Pkg	Price
GBD-S1	Glove Box Dispenser – Single	EA	\$53.00
GBD-T3	Glove Box Dispenser – Triple with Dividers	EA	\$95.00

More information available upon request.