


Date: 11/5/2018

Healthmark Industries introduces the Instrument Tube Brushes to its ProSys™ instrument care product line

Instrument Tube Brushes



Healthmark Industries is excited to announce the addition of the Instrument Tube Brushes to its ProSys™ instrument care product line. Designed as single-use brushes, the Instrument Tube Brushes assist with the cleaning of internal lumens and channels of various sizes. Manufactured to facilitate maximum cleaning performance while ensuring compatibility with the channels of cannulated items, the nylon bristle brush tips connect to a stainless-steel wire and joins to a versatile acetal tubing, which provides adequate debris removal without causing scope damage. Offered in lengths of thirty and sixty centimeters, the Instrument Tube Brushes have the following different colored shafts that indicate the brush head diameters: White (2mm & 15mm), Black (4mm), Blue (6mm), Green (9mm), and Gray (12mm). The ITB Brushes are available for purchase in a package of thirty.



healthmark
INDUSTRIES CO.
health care products

33671 Doreka
Fraser, Michigan 48026
Telephone: (586) 774-7600
Fax: (586) 774-6473
www.HMARK.com

Part#	Description	Pkg	Price
ITB-30-2	Cleaning Brush 2mm Dia., 30cm L	PK-30	\$57.75
ITB-30-4	Cleaning Brush 4mm Dia., 30cm L	PK-30	\$25.20
ITB-30-6	Cleaning Brush 6mm Dia., 30cm L	PK-30	\$31.50
ITB-30-9	Cleaning Brush 9mm Dia., 30cm L	PK-30	\$39.38
ITB-30-12	Cleaning Brush 12mm Dia., 30cm L	PK-30	\$49.22
ITB-30-15	Cleaning Brush 15mm Dia., 30cm L	PK-30	\$61.62
ITB-60-2	Cleaning Brush 2mm Dia., 60cm L	PK-30	\$64.05
ITB-60-4	Cleaning Brush 4mm Dia., 60cm L	PK-30	\$31.50
ITB-60-6	Cleaning Brush 6mm Dia., 60cm L	PK-30	\$44.10
ITB-60-9	Cleaning Brush 9mm Dia., 60cm L	PK-30	\$55.13
ITB-60-12	Cleaning Brush 12mm Dia., 60cm L	PK-30	\$68.91
ITB-60-15	Cleaning Brush 15mm Dia., 60cm L	PK-30	\$86.13

More information available upon request.