

Date: 6/26/2018

Healthmark Industries introduces the Loupe Set with Smartphone Adapter to its ProSys™ Optical Inspection product line

Loupe Set with Smartphone Adapter



Healthmark Industries is delighted to introduce the Loupe Set with Smartphone Adapter to its ProSys™ Optical Inspection product line. Manufactured for close visual inspection of medical devices, the ergonomically designed 4-piece eye loupe set and smartphone adapter are made from ABS plastic and rubber. Engineered to rest comfortably in your eye socket, the eye loupes can be used manually by hand or with a phone by connecting the eye loupes to the smartphone adapter and attaching it to the cellular device by placing it in front of the camera to capture pictures and videos for documentation purposes. The unique four eye loupes come in the following levels of magnification: 4.5x, 6.5x, 8x and 13x. The Loupe Set with Smartphone Adapter comes as complete set, with the 4 loupes and the Smartphone Adapter.

Part#	Description	Pkg	Price
LUP-004	Loupe Set with Smartphone Adapter	EA	\$30.00

More information available upon request.