

Date: 6/1/18

Healthmark Industries introduces the ChannelCheck™ Convenience Pack to its Cleaning Verification Product Line

ChannelCheck™ Convenience Pack



Healthmark Industries is delighted to introduce the ChannelCheck™ Convenience Pack to its Cleaning Verification product line. Designed as a simple flush method test, the ChannelCheck™ allows your department to simultaneously check for protein, carbohydrates, and hemoglobin residue after cleaning. The ChannelCheck™ is the first product that is capable of testing virtually any lumened instrument for residual organic soils, no matter the channel size. ChannelCheck™ helps ensure proper cleaning and helps to reduce potential risk to personnel and patients. The convenience pack consists of the following: 50 UCC-101 Test Strips (individually packaged), 1 ATS Control, 50 Zip Lock Bags for collection, 50 Prefilled 10ml Sterile Water Syringes, 1 UCC Color Chart, and 1 IFU.

Part#	Description	Pkg	Price
UCC-101CP	ChannelCheck Convenience Pack	EA	\$175.00

More information available upon request.