


|  |  |
|--|--|
| <b>Brand Name of Product</b>             | Long Cynch-Loks  |
| <b>Generic Name of Product</b>           | Tamper Evident Seals   |
| <b>Product Code Number(s)</b>            | 927-15, 927-07, 927-11, 1515 GN HLD  |
| <b>Intended Use:</b>                     | To show tampering and/or access to critical equipment and supplies.  |
| <b>Range of Applications for Product</b> | Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.   |
| <b>Key Specifications of Product</b>     | <ul style="list-style-type: none"> <li>• Numbered</li> <li>• Seal and broken by hand</li> <li>• Cannot be resealed once broken</li> <li>• Variety of colors available</li> </ul> |

**Shipping & Storage**

|   |  |
|---|--|
| <b>Shipping Conditions &amp; Requirements</b> |  |
| <b>Storage Conditions</b>                     |  |
| <b>Packaging Contents</b>                     |  |
| <b>Shelf Life</b>                             |  |

**Instructions for Using Product**

|  |  |
|--|--|
| <b>Description of Use (s)</b>            | Seal crash carts, emergency boxes and other areas where evidence of tampering is critical  |
| <b>Preparation for Use</b>               |  |
| <b>Diagrams (drawings, pictures)</b>     |   |
| <b>Steps for Use of Product</b>          | <ol style="list-style-type: none"> <li>1. Place strap through locking chamber with teeth towards neck of bag.</li> <li>2. Pull strap manually or place hook of puller tool in slotted strap-end. Pull tightly.</li> <li>3. For tight grip, place jaws of the tool over teeth &amp; pull. The bag is now securely sealed. To remove surplus strap, place jaws of tool over teeth at locking chamber &amp; twist strap.</li> </ol> |
| <b>Interpretation of Results</b>         |  |
| <b>Contraindications of Test Results</b> |  |
| <b>Documentation</b>                     | Follow the policy of your facility for recording date, staff and other information critical to management of critical inventory.   |
| <b>Special Warnings and Cautions</b>     | If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory.  |
| <b>Disposal</b>                          |  |

**Reprocessing Instructions**

|   |  |
|---|--|
| <b>Point of Use</b>                         |  |
| <b>Preparation for Decontamination</b>      |  |
| <b>Disassembly Instructions</b>             |  |
| <b>Cleaning – Manual</b>                    |  |
| <b>Cleaning – Automated</b>                 |  |
| <b>Disinfection</b>                         |  |
| <b>Drying</b>                               |  |
| <b>Maintenance, Inspection, and Testing</b> |  |
| <b>Reassembly Instructions</b>              |  |
| <b>Packaging</b>                            |  |
| <b>Sterilization</b>                        | Compatible with steam sterilization. Standard cycle: 132°C/4 minutes |
| <b>Storage</b>                              |  |
| <b>Additional Information</b>               |  |

|  |   |
|--|---|
| <b>Related Healthmark Products</b>     |   |
| <b>Other Product Support Documents</b> | Security Products Brochure, Security Products Price List  |
| <b>Reference Documents</b>             |   |
| <b>Customer Service Contact</b>        | Healthmark Industries Company, Inc<br>18600 Malyn Blvd.<br>Fraser, MI 48026<br>1-586-774-7600<br><a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a><br><a href="http://hmark.com">hmark.com</a> |

2020-10-08 msmith