



**Instructions for Use: Pouch Sterilization Rack**

<b>Brand Name of Product</b>	Pouch Sterilization Rack
<b>Generic Name of Product</b>	Sterilization Rack
<b>Product Code Number(s)</b>	SR-048, SR-010, SR-040, SR-004, SR-025, SR-005, SR-008, SR-016, SR-020
<b>Purpose of Product</b>	Keeps pouches on edge preventing over stacking, facilitates sterilant contact and reduces drying time. All in accordance with AAMI ST79
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	Constructed of stainless-steel rod and tungsten inert gas welded for extra strength. Styles include: <ul style="list-style-type: none"> <li>• Full Loop</li> <li>• Peg</li> <li>• Multi-Height Racks</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	SRS's peel pouch racks are designed for optimal airflow for the sterilization cycle and allow for complete drying of peel pouches. We offer an entire line of peel pack sizes and configurations customized for your requirements.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures):</b>	N/A
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Place sterilization rack in the designated area of the sterilizer.</li> <li>2. Place each sealed pouch in each slot so there is space between each pouch.</li> <li>3. When sterilization is complete, remove each pouch from the rack.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	Wash with a mild detergent. Rinse thoroughly.
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	Compatible with steam, EO and gas plasma/VHP sterilization.
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	Sterilization Brochure, Sterilization Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026

	1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com
--	---