



Instructions for Use: Silicone Tray Corner Protector

Brand Name of Product	Silicone Tray Corner Protectors
Generic Name of Product	Tray Corner Protector
Product Code Number(s)	CP-80NS, 50062, 50062YL
Purpose of Product	To protect the corner of the sterilization wrap from the sharp edges of sterilization trays.
Range of Applications for Product	Wrapped trays that are steam, EO or gas plasma sterilized. These tray corners are reusable.
Key Specifications of Product	Material made of silicone. <ul style="list-style-type: none"> • CP-80NS: Pink • 50062: Purple • 50062YL: Yellow

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use (s)	These reusable tray corners are a perfect solution for problematic corner tears.
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	Place tray on the tray corners. Be sure that all four corners are covered.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Tray Corner Protector can be cleaned manually with detergent and water.
Cleaning – Automated	Tray Corner Protectors can be processed through an automated instrument washer.
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Compatible with steam, VH2O2/O3/EO sterilization.

Storage	N/A
Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	Sterilization Brochure, Sterilization Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2021-06-14 Matt Smith