

| | |
|--|---|
| Brand Name of Product | Sterimarker |
| Generic Name of Product | Sterimarker |
| Product Code Number(s) | MAR1, MAR2, MAR3, 1753DZ |
| Intent of Use | To write on packages, in accordance with the IFU of the packaging material. |
| Range of Applications for Product | The special ink is immediately dry, waterproof, lightproof, resistant to atmospheric action, heat resistant, and acid resistant. No chemical or toxic substances are released during use. The ink marking withstands sterilization by steam, EO, VH2O2 and radiation. |
| Key Specifications of Product | <ul style="list-style-type: none"> ● Nontoxic ● Waterproof ● No chemicals or toxic substances are released during use. |

| Shipping & Storage | |
|---|---|
| Shipping Conditions & Requirements | N/A |
| Storage Conditions | <ul style="list-style-type: none"> ● No special precautions are necessary regarding storage/transportation. ● Do not expose them to very high humidity or extreme heat. |
| Packaging Contents | N/A |
| Shelf Life | N/A |

| Instructions for Using Product | |
|--|---|
| Description of Use(s) | To write on packages in accordance with the IFU of the packaging material. |
| Preparation for Use | N/A |
| Diagrams (drawings, pictures) | N/A |
| Steps for Use of Product | <ol style="list-style-type: none"> 1. Uncap marker. 2. Write on surfaces identified as appropriate to write on, according to the product's instructions. 3. Recap marker after use. |
| Interpretation of Results | N/A |
| Contraindications of Test Results | N/A |
| Documentation | N/A |
| Special Warnings and Cautions | <ul style="list-style-type: none"> ● Recap marker immediately after use. ● The marker pen is not autoclavable. ● Do not write on surfaces with a glossy finish, the ink will smudge after sterilization. |
| Disposal | N/A |

| Reprocessing Instructions | |
|---|--|
| Point of Use | N/A |
| Preparation for Decontamination | N/A |
| Disassembly Instructions | N/A |
| Cleaning – Manual | N/A |
| Cleaning – Automated | N/A |
| Disinfection | N/A |
| Drying | N/A |
| Maintenance, Inspection, and Testing | N/A |
| Reassembly Instructions | N/A |
| Packaging | N/A |
| Sterilization | Ink is compatible with Steam, VH2O2/O3/EO sterilization and radiation. |
| Storage | N/A |
| Additional Information | N/A |

| | |
|------------------------------------|-----|
| Related Healthmark Products | N/A |
|------------------------------------|-----|

| | |
|--|--|
| Other Product Support Documents | N/A |
| Reference Documents | Sterilization Brochure, Sterilization Price List |
| Customer Service Contact | Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com |