



Instructions for Use: Count Sheet Holders

Brand Name of Product	Count Sheet Holders
Generic Name of Product	Count Sheet Holders
Product Code Number(s)	PBR-46, PBR-76
Purpose of Product	Designed to hold instrument count sheets on the outside of wrapped trays (in compliance with AORN recommended practices).
Range of Applications for Product	For all trays with count sheets that accompany the set of surgical instruments.
Key Specifications of Product	<ul style="list-style-type: none"> • PBR-46: <ul style="list-style-type: none"> ○ Transparent poly front with medical grade paper back. ○ 4- x 6 inches. • PBR-76: <ul style="list-style-type: none"> ○ Medical grade paper back ○ 7 ½- x 3 inches.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Holds instrument count sheets on the outside of wrapped trays.
Preparation for Use	N/A
Diagrams (drawings, pictures):	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Fold the count sheet over to fit inside the medical grade paper pouch. 2. Peel away the protective layer to expose adhesive backing and adhere the medical grade paper pouch to the side of your wrapped package.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	Compliant with AORN guidelines, which recommend placing count sheets outside of sterilization trays.
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Materially compatible with steam and EO sterilization.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	Slotted Tray Labels
Other Product Support Documents	Sterilization Packaging Price List, Sterilization Packaging Brochure

Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com