



**Instructions for Use: Count Sheet Holders**

<b>Brand Name of Product</b>	Count Sheet Holders
<b>Generic Name of Product</b>	Count Sheet Holders
<b>Product Code Number(s)</b>	PBR-46, PBR-76
<b>Purpose of Product</b>	Designed to hold instrument count sheets on the outside of wrapped trays, in compliance with AORN recommended practices.
<b>Range of Applications for Product</b>	All trays where a count sheet accompanies the set of surgical instruments.
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• PBR-46: <ul style="list-style-type: none"> <li>○ Medical grade paper</li> <li>○ 6" x 4"</li> </ul> </li> <li>• PBR-76: <ul style="list-style-type: none"> <li>○ Transparent poly front with medical grade paper back</li> <li>○ 7.5" x 3"</li> </ul> </li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	Come packaged in quantity of 1000.
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	Holds instrument count sheets on the outside of wrapped trays.
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ul style="list-style-type: none"> <li>• Simply fold the count sheet over to fit inside medical grade paper pouch.</li> <li>• Adhere to the side of your wrapped package.</li> </ul>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	Compliant with AORN guidelines which recommend placing count sheets outside of sterilization trays.
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	Materially compatible with steam and EO sterilization.
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	Slotted Tray Labels
<b>Other Product Support Documents</b>	Sterilization Packaging Price List, Sterilization Packaging Brochure
<b>Reference Documents</b>	
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600

	<a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com
--	---

2012-09-11 msmith