



healthmark
INDUSTRIES CO.
health care products

33671 Doreka
Fraser, Michigan 48026
Telephone: (586) 774-7600
Fax: (586) 774-6473
www.HMARK.com

Product Evaluation Form

Product: Steriking See Through Sterilization Pouches:

Steriking pouches provide optimum protection and microbial barrier for sterile contents during storage, handling, transportation and aseptic presentation. Steriking pouches are constructed of high weight medical grade paper with a five layer multi film, making these pouches the most durable and toughest available. Product range provides a wide variety of sizes and types for optimum choice and correct pack sizing. Pouches are available on rolls, flat or gusseted, individual sized, with heat or self-seal options to meet individual user needs. Pouches are validated for double pouching when appropriate.

Evaluator's Name: _____ Date _____
 Title: _____ Dept: _____
 Phone: _____

Hospital/Facility: _____ Healthmark Rep: _____
 Address: _____ City _____ State: _____

Evaluation Start Date: _____ Evaluation End Date: _____

Product type circle which applies: *heat seal self seal tyvek roll gusseted*

Using the scale below, please rate the following features:

4 = Excellent 3 = Good 2 = Fair 1 = unacceptable NA = Non applicable

- Heat Seal Pouch* - film turns to darker green allows visual check of heat seal integrity _____
- Self Seal Pouch* - Strong secure seal closure _____
- Clean peel fiber free opening for sterile presentation _____
- Strong seal closure and seams prevents bursting _____
- Sealed edges in upper corners prevents dust reducing risk of contamination when opening _____
- Ease of use and handling -thumb notches at pouch ends facilitates ease in filling & opening _____
- Safe print design indicators and imprints on outside avoids risk of ink contamination _____
- Traceability each pouch contains a manufacturers lot code _____
- Availability and value of manufacturer's product information and validation _____
- Process indicators produce clearly contrasting color change to distinguish processed goods _____
- Multi layer film and high weight paper provide superior barrier, durability prevents tears and punctures _____
- Types, sizes and options meet varied and unique needs _____
- I would recommend the use of this product _____

Comments: _____