

## Change in STERIKING® TYVEK® packages design

The STERIKING® LTR Rolls are having the symbol meaning “Do not use if sterile barrier system is damaged”. Now the symbol (in Figure 1) is taken in to use also in STERIKING® Tyvek® pouches (LTS- and LTSS-pouches) where it replaces the text in different languages.




Figure 1

In Tyvek® pouches the logo VH2O2 replaces the phrase “processed in plasma” like it has been in Tyvek® rolls. This logo indicates that STERIKING® Tyvek® products can be used in other hydrogen peroxide sterilization processes as well than gas plasma.

At the same time into the design of both Tyvek® pouches and rolls is added additional text “DuPont Tyvek®” to clarify that the Tyvek® used also in STERIKING® Tyvek® products is made by DuPont.

Changes are valid for pouches manufactured from December 2011 onwards.



\*) Note: The colored frames in this illustration highlight the text blocks for further clarification and are not considered as a part of the design

Sincerely,  
**WIPAK OY**