



SPECIFICATION

REF 01413	Product: Indicator dots with indicator	Size 10,0 x 10,0 mm
-----------	---	------------------------


Material: Tyvek self-adhesive,

Indicator: H₂O₂ (Plasma) indicator yellow with H₂O₂ indicator according to DIN EN ISO 11140-1 class colour change from yellow / light orange to blue


Fabrication: 25 pcs./sheet
label size: 10 x 10 mm

Packaging: 20 sheets = 500 pcs. per PE-bag


Label:



Layout:
before sterilization



after sterilization



Distributed by



healthmark
 Industries Co.
 health care products

33671 Doreka Drive Fraser, MI 48026
 1-800-521-6224 www.hmark.com