



## Product Technical Data Sheet

Form No.: FTP-094 Rev.: B

**Product Name:** SteriTec EO Indicator Labels

**Product Type:** Process Indicator

**SteriTec Reorder Number:** CI 126 Rev. A

**Description:** The EO Indicator labels are designed for use as process indicators for industrial ethylene oxide gas sterilization processes.

### Performance Characteristics:

The indicator label will change from yellow to blue when exposed to EO gas sterilization conditions.

**Color Change:** yellow to blue

**Complete Color Change Time:** 15 minutes @ 130° F, 600mg/l EO gas, 40-50% relative humidity.

**Indicator Stability:** Excellent before and after sterilization

### Physical Characteristics:

**Background Print Color:** N/A

**Dimensions:** ½" diameter

**Substrate Material:** Label Stock

**Units / Case:** 10 rolls

**Selling Unit:** Case

### Environmental Limitations:

**Storage:** 50°-100° F (10°-38° C) at 30-60% humidity.

**Shelf Life:** 3 years

**SteriTec FDA 510(k) Release to Market Date:** None. Industrial Use only.

### Standards Classifications:

EN 867-1, Class A

ISO 1140-1, Class 1