



## Product Technical Data Sheet

Form No.: FTP-170 Rev.: D

### Product Name: SteriTec Universal Container Card

**Product Type:** Process Indicator

**SteriTec Reorder Number:** CI 138

**Description:** The SteriTec Universal Container Card is a steam sterilization process indicator designed for use with various types of sterilization containers. The Container Cards will function in sterilizers that operate at 121°C or 134°C. It provides a visible confirmation that the instrument container being sterilized has been exposed to steam sterilization processing.

#### Performance Characteristics:

When exposed to steam sterilization conditions, the bar changes color.

**Color Change:** White to Black

**Complete Color Change Time:** 10 minutes +/- 30 seconds at 121°C  
2 minutes +/- 30 seconds at 132°C-134°C

**Indicator Stability:** Excellent Before and After Exposure

#### Physical Characteristics:

The card can fit in various container card-holders by folding or not folding the top sizing flap to the appropriate size.

**Background Print Color:** Black/Gray

**Dimensions:** 1.531" W X 2.75" L or with flap folded down 1.4" W X 2.75" L

**Substrate Material:** 80 lb Cover

**Units / Box:** 500

**Boxes/Case:** 6

**Selling Unit:** Case

#### Environmental Limitations:

**Storage:** 50°–100°F (10°–38°C) at 10-70% relative humidity

**Shelf Life:** The expiry date is printed on the product packaging.

**SteriTec FDA 510(k) Release to Market Date:** K052516.

#### Standards Classifications:

EN ISO 11140-1, Class 1

DISTRIBUTED BY:  
HEALTHMARK INDUSTRIES, CO.  
33671 DOREKA DRIVE  
FRASER, MI 48026  
1-800-521-6224 / WWW.HMARK.COM