



Healthmark Industries Company

BI Interpretation Guide

INSTANT READ-OUT BIOLOGICAL

Steam Sterilization Monitor



FAIL

PASS

SteriTec INSTANT READOUT INTEGRATOR
 STERILIZER NO.: _____ LOAD NO.: _____
 DATE: _____

PASS 270 °F
 IF GREEN (COLOR STANDARD)

RESULTS: STERILIZATION CONDITIONS
 ACHIEVED DATE: _____
 NOT ACHIEVED RELEASED BY: _____

SteriTec

SteriTec INSTANT READOUT INTEGRATOR
 STERILIZER NO.: _____ LOAD NO.: _____
 DATE: _____

PASS 270 °F
 IF GREEN (COLOR STANDARD)

RESULTS: STERILIZATION CONDITIONS
 ACHIEVED DATE: _____
 NOT ACHIEVED RELEASED BY: _____

SteriTec

INSTANT READ-OUT BIOLOGICAL

Ethylene Oxide Sterilization Monitor



FAIL

PASS

INSTANT READOUT **EO** INTEGRATOR
 STERILIZER NO.: _____ LOAD NO.: _____
 DATE: _____

PASS (COLOR STANDARD)

RESULTS: STERILIZATION CONDITIONS
 ACHIEVED DATE: _____
 NOT ACHIEVED RELEASED BY: _____

SteriTec

INSTANT READOUT **EO** INTEGRATOR
 STERILIZER NO.: _____ LOAD NO.: _____
 DATE: _____

PASS (COLOR STANDARD)

RESULTS: STERILIZATION CONDITIONS
 ACHIEVED DATE: _____
 NOT ACHIEVED RELEASED BY: _____

SteriTec