



healthmark
INDUSTRIES CO.
health care products

33671 Doreka
Fraser, Michigan 48026
Telephone: (586) 774-7600
Fax: (586) 774-6473
www.HMARK.com

UNDERGUARD LINERS (SOAKER SHEETS)

Made of the same pure cellulose as our OrthoShield™ sheets, these UnderGuard™ sheets are sized for use on autoclave loading carts.

These sheets serve to absorb and rapidly disperse condensate during sterilization. Further, they will help prevent dripping of condensate between shelves.

Best of all, they protect wrapped trays for any sharp edges that might tear wrap when trays are loaded or unloaded.

Liners come in three standard sizes to fit both floor loaded and shelf loading autoclaves.

Code Number	Width	Length	Case Quantity
SOK008SS	22.00	48.00	100
SOK009SS	22.00	60.00	100
SOK010SS	20.00	36.00	50

PLEASE CALL TODAY FOR YOUR PRICING

1-800-521-6224

WWW.HMARK.COM