


<b>Brand Name of Product</b>	UnderGuard for Plasma™
<b>Generic Name of Product</b>	Tray liners
<b>Product Code Number(s)</b>	MT-2030, MT-1520, MT-1224.
<b>Purpose of Product</b>	For use under the tray to protect outer wrap from sharp corners and imperfections of sterilization trays.
<b>Range of Applications for Product</b>	
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Manufactured from spun polypropylene.</li> <li>The unique purple color helps distinguish this protective layer from sterilization wrap.</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	+15°C-30°C and 30-70% RH
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use (s)</b>	To protect outer wrap from sharp corners and imperfections of sterilization trays.
<b>Preparation for Use</b>	<ol style="list-style-type: none"> <li>Select a liner which will fit best the target basket.</li> <li>If using with wrapped tray, liner must be inside the wrap and may be used inside or under the basket.</li> <li>If using with a closed container, liner must be used inside the container, either inside or below the basket.</li> </ol>
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>Place liner either inside or beneath basket. In the event you are using a container, the liner should be placed in case first.</li> <li>Place items in tray.</li> <li>Use proper sterility assurance indicators and follow hospital protocol.</li> <li>Wrap or close package.</li> <li>Place in gas-plasma/VHP or steam sterilizer.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	MATxxx are number 5 for recycling.

Reprocessing Instructions	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	

<b>Packaging</b>	
<b>Sterilization</b>	Material compatible with gas-plasma/VHP and steam sterilization up to 137°C
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Sterilization Brochure, Sterilization Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2018-07-09 msmith